

Upper Anteriors

Upper Anteriors

Lower Anteriors


Upper and Lower Posteriors

Choice System

Choice in relation to men size

E1

A 43,3
B 8,7
C 8,3




F4

A 48,3
B 10,4
C 8,8




E2

A 42,7
B 9,2
C 8,1




G1

A 40,8
B 8,0
C 7,6




E3

A 45,0
B 10,1
C 9,0




G2

A 44,5
B 9,2
C 8,6




E4

A 48,8
B 10,5
C 9,0




G3

A 44,8
B 9,8
C 8,5




H1

A 40,6
B 8,9
C 7,4




H4

A 43,3
B 9,0
C 8,0




K1

A 32,3
B 8,1
C 4,4



K3

A 33,8
B 8,5
C 4,5



K6

A 35,5
B 9,8
C 5,2



K9

A 38,5
B 9,5
C 5,6



M1

A 59,2
B 61,7



M2

A 59,4
B 62,8



M4

A 62,6
B 67,3




M5


A 67,0
B 71,5






leptosomal



pyknic




athletic




Anterior Moulds Classification

Shade Guide and Combination Chart


Structure



Mesio-Distal Profile




curved




flat

Lateral Profile



curved



flat

A1	B1	C1	D2
A2	B2	C2	D3
A3	B3		
A3,5			

Uppers	Lowers	Posteriors
E1	K3	M2
E2	K6	M2
E3	K6	M4
E4	K9	M5
F4	K9	M5
G1	K1	M1
G2	K6	M2
G3	K6	M4
H1	K1	M1
H4	K3	M4

Craniform

Cranial GmbH

Cranial GmbH
 Telefon: +49 (0)208 - 488 37 21 · Fax: +49 (0)208 - 488 37 23
 Mail: info@cranial.de · www.cranial.de

by **Cranial** GmbH